

REF LS 2000-12 / LS 2000-13

General Description

FLEX LITE series provides deep cavity lighting, direct illumination into the surgical field. Achieves optimal visual assistance, eliminating shadows.

It has a bendable and malleable strip with 3M double coated medical tape that can be easily shaped and attached to most retractors or instruments.



Read the operating instruction carefully before use. Failure to follow may cause damage or injuries.



Warnings

Failure to follow may cause harm to person.

- This product is not recommended for direct contact with the central nervous system.
- Avoid getting in contact with massive fluid for prolonged time.
- Sterile only if package is unopened and undamaged.
- DO NOT operate damaged or malfunctioning instrument.
- Single use only! DO NOT Re-sterilize!
-

Caution

Failure to follow may cause damage to product.

- DO NOT drop or hit the instrument.
- Store in a cool dry place.
- No modification of this equipment is allowed.
- Do not mix with used or other types of batteries.

Indications


Provide visible illumination of the surgical field or the patient.

Contraindications

None


Instructions for Use

1. Peel open the pack and remove all of the items with aseptic technique.
2. Tear to release liner from 3M double coated medical tape.
3. Bend and attach Light strip (A) on the retractor.
4. Depress the button (B) to switch ON for use.
5. Depress the button (B) again to switch OFF.



Battery Removal

1. Use a tool (e.g. flat screwdriver) to open from the slot behind the housing.
2. Cut off the positive and negative terminals.



3. Place a masking tape on both positive and negative terminal ends to ensure the batteries are properly isolated from metal or other batteries to avoid a short circuit.

Disposal

- The empty batteries should be disposed of through specially designated collection boxes, recycling points or electronics retailers.
- DO NOT dispose battery in fire or recharge.
- Dispose the device as appropriate, immediately after use.
- Ensure that the device is disposed of safely and in accordance with all current national and international waste disposal directives.

Shelf Life

Please check package label.

Transportation and Storage

- Temperature: -25~70°C
- Relative Humidity: 10%~95%RH maximum
- Atmospheric Pressure: 860hPa~1060hPa


















Operation

- Operating temperature: 5~40°C
- Relative Humidity: 15%~85%RH maximum
- Atmospheric Pressure : 860hPa~1060hPa

Specification

- AAA 1.5V x 2 Alkaline Batteries.
- Length of Light Strip: 100 mm (LS2000-13) or 200mm (LS2000-12).
- Typical color temperature 5700~6000K.
- 3M™ Double Adhesive Medical Tape.
- The device complies with the following standards EN 60601-1-2, EN 60601-1, EN 62471.

Labeling Symbols

	Refer to instruction manual/booklet
	General warning sign / Caution
	Do not use if package is damaged
	Do not re-use
	Keep dry
	Temperature limit
	Sterilized using ethylene oxide
	Catalog number
	Batch code
	Serial number
	Complies with EU directives
	Date of manufacture
	Manufacturer
	Authorized representative in the European Community
	Use by date
	Liquid ingress protection: immersion up to 1m depth
	Recycle alkaline battery

